Berichrom* C1-Inhibitor

Reagents for the determination of C1 inhibitor

Intended Use

Determination of the functional activity of C1-inhibitor in plasma for the diagnosis of diminished C1-inhibitor synthesis, increased consumption and for monitoring substitution therapy and steroid therapy.



Summary and Explanation

C1-inhibitor is the only inhibitor for the activated complement components C1r and C1s which occur at the start of the classical pathway of complement activation. C1-inhibitor therefore prevents excessive activation of the complement system. In addition, C1-inhibitor is physiologically the most important plasmatic inhibitor for activated coagulation factor XIIa (FXIIa) and for kallikrein, rapidly and irreversibly forming an inactive complex with each. Thus, C1-inhibitor is the most important regulator of the contact phase¹.

Hereditary deficiency of C1-inhibitor causes hereditary angioneurotic oedema (HANE), a special form of Quincke's oedema². Acquired deficiency states occur in malignant diseases, liver cirrhosis, preeclampsia, pneumonia and other infectious diseases. The determination of C1-inhibitor is furthermore indicated for monitoring substitution therapy with C1-inhibitor as well as for monitoring steroid therapy.

Principle of the Method

C1-inhibitor in the sample inhibits a receiving volume of C1-esterase. The remaining activity of the C1-esterase is determined by measuring the increase in absorbance at 405 nm in a kinetic test based on the following reactions:

Reagents

Materials provided

Berichrom* C1-inhibitor, Testkit, Code No. OUIA with

3 x for 5 mL C1-Esterase

3 x for 1 mL Substrate Reagent

1 x 3 mL Substrate Solvent

Composition

C1-Esterase (human), lyophilized; concentration in the working solution:

KH₂PO₄ 65 mmol/L, lysin 1 g/L, mannitol 5 g/L, pH 7.5

Preservatives:

5-chloro-2-methyl-4-isothiazole-3-one (max. 8 mg/L) and 2-methyl-4-iso-thiazole-3-one (max. 2 mg/L)

Substrate Reagent: lyophilized; concentration in the working solution:

methylloxycarbonyl-L-lysyl(ε-carbobenzoxyl)-Glycyl-L-Arginyl-p-nitrolanilide (MEOC-Lys(ε-Cbo)-Gly-Arg-pNA) 5 mmol/L, mannitol 20 g/L.

Substrate Solvent: Polypropylene glycol, 100 g/L

Warnings and Precautions

1. For *in vitro* diagnostic use.

2. As complete absence of infectious agents can never be assured, all materials derived from human tissue or body fluids should always be handled with due care, observing the precautions recommended for biohazardous material³.

Preparation of the Reagent

C1-Esterase (human): Dissolve vial contents in the labelled amount of distilled water.

Substrate Reagent: Dissolve vial contents in the labelled amount of substrate solvent.

Mix carefully once more before using.

Storage and Stability

Store the test kit unopened at +2 to +8 °C and use by the date given on the label.

Stability after reconstitution:

Temperature	C1-Esterase	Substrate Reagent
+37 °C	8 hours	1 week
+30 °C	1 day	2 weeks
+15 to +25 °C	1 week	2 weeks
+2 to +8 °C	2 weeks	6 weeks
-20 °C	1 month	6 months

Stored in their original vials the C1-Esterase can be frozen up to 3 times, the substrate up to 10 times. Different coagulation analyzers have individual stability data.

Materials required but not provided

Standard Human Plasma, Code No. ORKL

Control Plasma N, Code No. ORKE

Control Plasma P, Code No. OUPZ

Acetic acid 20% (only for the two-point method)

Imidazole buffer solution, Code No. OQAA or

Dade® Owren's Veronal Buffer, Code No. B4234 or

Dade® CA System Buffer, Code No. B4265 or

Physiological saline solution

Coagulation analyzer (see Chapter titled "Limitations of the Procedure")

Equipment

Berichrom* C1-Inhibitor can be used in a great number of coagulation analyzers.

Observe the operating instructions issued by the analyzer's manufacturer!

Specimens

To obtain the plasma, carefully mix 1 part sodium citrate solution (0.11 mol/L) with 9 parts venous blood, avoiding the formation of foam. Centrifuge immediately at no less than 1500 x g for at least 10 minutes.

Stability of the sample:

-20 °C 1 month

+ 2 to + 8 °C 2 days

+15 to +25 °C 6 hours

Plasma stored at -20 °C is to be thawed within 10 minutes at 37 °C after which the assay is to be performed within 2 hours. In some cases the assay can be performed using serum, however the plasma and serum values do not always correlate.

Procedure

Manual:

Semi-micro test:

cuvette: 1 cm path length

wavelength: 405 nm

test temperature: optionally +25 °C, +30 °C or +37 °C

Pre-warm the C1-esterase, substrate reagent and plastic cuvettes/tubes to the selected test temperature.

Pipetting scheme for the kinetic method

	C1-esterase enzyme value (CEV)	Sample
Isotonic saline	20 µL	–
Plasma sample	–	20 µL
C1-esterase	1000 µL	1000 µL
Mix and incubate at selected test temperature for exactly 5 min.		
Substrate Reagent	100 µL	100 µL
Mix and determine ΔA _{405 nm} /minute		
Photometer/Stopwatch: Read the absorbance within 30 seconds and start the stopwatch simultaneously. Read the absorbance again after exactly 60 and 120 seconds, and calculate the respective ΔA/minute values by subtraction, and calculate the mean of the two values.		

Test mixtures for the Chromotimer: Use the same incubation times but half the volumes of sample and reagents.

Evaluation, kinetic method

Each assay series requires at least one C1-esterase enzyme value (ΔA/min_{CEV}) and one reference measurement value, which is determined using a reference plasma with an assigned value for C1-inhibitor (e.g. Standard Human Plasma) as the sample. These values are used to calculate the laboratory-internal factor F_L:

$$F_L = \frac{\text{Assigned value (\% of the norm)}_{\text{reference plasma}}}{\Delta A / \text{minute}_{\text{CEV}} - \Delta A / \text{minute}_{\text{reference plasma}}}$$

The C1-inhibitor content of the sample in % of the norm can then be calculated from the following formula:

$$\text{C1-inhibitor}_{\text{sample}} (\% \text{ of the norm}) = F_L \times (\Delta A / \text{minute}_{\text{CEV}} - \Delta A / \text{minute}_{\text{sample}})$$

Pipetting scheme for the two-point method

Pre-warm the C1-esterase, substrate reagent and plastic cuvettes/tubes to +37 °C.

	C1-esterase enzyme value (CEV)	Sample	Sample blank
Isotonic saline	20 µL	–	1100 µL
Plasma sample	–	20 µL	20 µL
C1-esterase	1000 µL	1000 µL	–
Mix and incubate for exactly 5 min at +37 °C.			
Substrate Reagent	100 µL	100 µL	–
Mix immediately and start the stopwatch simultaneously. After exactly 2 min, add:			
Acetic acid 20%	500 µL	500 µL	500 µL
Mix immediately and read the absorbance against the sample blank within 60 min.			

Evaluation, two-point method

Each assay series requires at least one C1-esterase enzyme value (A_{CEV}) and one reference measurement value, which is determined using a reference plasma with an assigned value for C1-inhibitor (e.g. Standard Human Plasma) as the sample. These values are used to calculate the laboratory-internal factor F_L:

$$F_L = \frac{\text{Assigned value (\% of the norm)}_{\text{reference plasma}}}{A_{\text{CEV}} - A_{\text{reference plasma}}}$$

The C1-inhibitor content of the sample in % of the norm can then be calculated from the following formula:

$$\text{C1-inhibitor}_{\text{sample}} (\% \text{ of the norm}) = F_L \times (A_{\text{CEV}} - A_{\text{sample}})$$

Notes

1. Doubling or halving the volumes used in the test mixtures has no effect on the calculation or laboratory-internal calculation factor.
2. F_L is a calculation factor established within the laboratory and needs to be determined only once in an assay series. This method of evaluation is not dependent on the selected test temperature and is basically equivalent to an evaluation via a reference curve.
3. If the sample has a C1-inhibitor content above 150% of the norm, the assay should be repeated using a sample prediluted 1+1 with isotonic saline. In this case the result is then multiplied by 2. If the sample has a C1-inhibitor content below 20% of the norm it is recommended to use 40 µL sample in the test. The result is divided by 2.
4. The laboratory-internal factor F_L or the calibration curve must be re-determined for each change of device and for each new lot of Berichrom* C1-Inhibitor.

Automatic:

Special instructions for coagulation analyzers may be obtained from Dade Behring by request.

Internal quality control

Normal range: Control Plasma N

Pathological range: Control Plasma P

Two controls should be measured with each calibration and at least every 8 hours during each testing day (one in the normal range and one in the pathological range). The controls should be processed just like the samples. Each laboratory should determine its own quality control range, either by means of the target values and ranges provided by the manufacturer of the controls or by means of its own confidence level established in the laboratory. If the measured control value lies outside the confidence level previously established, then the reagents, the laboratory-internal factor F_L and/or the calibration curve and the coagulation analyzer should be examined.

Limitations of the Procedure

There are appropriate coagulation analyzers for which Dade Behring provides operating instructions. The utilization of Berichrom* C1-Inhibitor in other coagulation analyzers must be validated by the testing laboratory. In doing so, performance data other than that provided may be obtained. Interference may be caused by high levels of lipids in the plasma samples.

Reference Interval

70–130%²

Systematic deviations from these ranges may be determined by the particular device. It may be necessary to calculate your own reference range in the laboratory.

Specific Performance Characteristics

Precision

The precision of Berichrom* C1-Inhibitor on the BCT* was estimated with Control Plasma N and Control Plasma P from a total of 8 separate runs over 5 testing days. The intra-assay coefficients of variation (CV) ranged between 1.8% and 7.9% and the inter-assay CV were between 3.2% and 6.6%.

Bibliography

1. Harpel P C: Blood Proteolytic Enzyme Inhibitors, in: Haemostasis and Thrombosis, Colman R W et al. (eds), Lippincott, Philadelphia (1987), 219–234
2. Thomas L: Komplement-System, in: Labor und Diagnose (4. erw. Auflage), Med. Verl.-Ges., Marburg (1992), 964–980
3. U.S. Department of Health and Human Services CDC, Biosafety in Microbiological and Biomedical Laboratories, HHS Publication (CDC) 93-8395; 1999; Section II; 8-16

Dade is a registered trademark of Dade Behring, Inc. in the USA, Germany and other countries.

* Berichrom and BCT are registered trademarks of Dade Behring Marburg GmbH in Germany and other countries.

 Dade Behring Marburg GmbH
Emil-von-Behring-Str. 76
D-35041 Marburg
www.dadebehring.com

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Dade ist eine eingetragene Marke der Dade Behring Inc. in den USA, Deutschland und anderen Ländern.
* Berichrom und BCT sind eingetragene Marken der Dade Behring Marburg GmbH in Deutschland und anderen Ländern.