N Antisera to Human Complement Factors (C3c, C4)

Intended Use

In vitro diagnostic reagents for the quantitative determination of complement factors (C3/C3c and C4/C4c) in human serum or heparinized or EDTA plasma by means of immunonephelometry on the BN[®] Systems as an aid in the diagnosis of immunologic disorders associated with complement C3 or C4 protein.

Summary and Explanation

The complement system is an integral part of the antigen-nonspecific immune defense. It can be activated via two reaction pathways, the classical pathway which is triggered primarily by cell-bound immune complexes, and the alternative pathway which is activated primarily by foreign bodies such as microorganisms. The complement component C3 is a key protein in both reaction pathways whereas C4 belongs to the classical pathway of complement activation. Complement activation is associated with consumption of components C3 or C4 so that a reduction in their concentrations can allow diagnostic conclusions to be reached. Diminished serum concentrations of C3 and C4 are observed primarily in active systemic lupus erythematosus (SLE), in forms of membrane proliferative glomerulonephritis and in immune complex diseases (serum sickness). In the case of SLE the serum concentrations of the complement factors reflect the activity of the disease. Diminished C3 values occur in acute glomerulonephritis and in membrane proliferative glomerulonephritis whereas isolated diminished levels of C4 can occur in hereditary angioedema (HAE) and in cases of autoimmune hemolytic anemia. Both complement components react as acute-phase proteins and may therefore show elevated serum concentrations in patients with inflammatory diseases. Hereditary deficiency states of both complement factors have been reported¹⁻³.

Principle of the Method

Proteins contained in human body fluids form immune complexes in an immunochemical reaction with specific antibodies. These complexes scatter a beam of light passed through the sample. The intensity of the scattered light is proportional to the concentration of the relevant protein in the sample. The result is evaluated by comparison with a standard of known concentration.

Reagents

Materials provided

N Antiserum to Human C3c, Code No. OSAP *or*
N Antiserum to Human C4, Code No. OSAO
one vial each, containing 5 mL or 2 mL

Composition

N Antisera are liquid animal sera and are produced by immunization of rabbits with highly purified human complement factors (C3c or C4). The antibody titers (T) are determined by radial immunodiffusion and are printed on the vial labels. The titers indicate the quantity of antigen in mg which will be precipitated in an agarose gel by 1 mL of the corresponding antiserum.

Preservative

N Antiserum to Human C3c or C4: sodium azide: <1 g/L

Warnings and Precautions

- For *in vitro* diagnostic use.
- Reagents containing sodium azide must be handled with due caution:
Do not ingest or allow to contact skin or mucous membranes! Sodium azide can form explosive azides when contacting heavy metals such as copper or lead.

Preparation of Reagents

The N Antisera are ready-for-use as supplied and require no additional preparation.

Storage and Stability

Stability at +2 to +8 °C:

See expiry date on label;

Stability once opened:

Four weeks if stored at +2 to +8 °C securely capped immediately after each use and contamination (e.g., by microorganisms) is precluded;

During storage, N Antisera can develop precipitates or turbidity which are not caused by microbial contamination and which do not affect their activity. In such cases, the antiserum should be filtered prior to use. Disposable filters with a pore size of 0.45 µm are suitable for this purpose.

Do not freeze.

On-board stability:

A minimum of five days at eight hours per day for 5 mL vials, and three days at eight hours per day for 2 mL vials or a comparable period of time.

Note: On-board stability may vary, depending on BN[®] System used and laboratory conditions. For further details, refer to BN[®] II or BN ProSpec[®] System Instruction Manual.

Materials required but not provided

BN[®] System

N Protein Standard SL (human), Code No. OQIM

N/T Protein Control SL/L, M and H (human), Code Nos. OQIN, OQIO and OQIP

N Reaction Buffer, Code No. OUMS

N Diluent, Code No. OUMT

BN[®] II Evaporation Stoppers (optional), Code No. OVLE

Additional materials and supplies as described in your BN[®] System Instruction Manual.

Specimens

Suitable samples are human serum or heparinized or EDTA plasma. Serum samples must be completely coagulated and, after centrifugation, must not contain any particles or traces of fibrin. Lipemic samples or frozen samples which became turbid after thawing must be clarified by centrifugation (10 minutes at approximately 15000 x g) prior to testing.

In the assay of C3, one should note that the antiserum is directed against the C3c fragment of the C3 molecule. The extent of fragmentation of C3 to the C3c fragment varies depending on the age of the sample and storage conditions. For fresh samples the C3 values obtained in the immunonephelometric assay have been observed to be up to 30 % lower than those obtained for stored samples, depending on the extent to which fragmentation has advanced⁴. Fragmentation of C4 to the more stable C4c fragment proceeds more slowly during storage of the sample than fragmentation of C3. Here, too, the C4 values obtained in the immunonephelometric assay have been observed to increase in relation to the period of storage. During storage for eight days at +2 to +8 °C or for three months at below -20 °C, serum or heparinized or EDTA plasma specimens may increase in C3c concentration by up to 17 %. C4c may increase under these storage conditions by up to 8 %. Therefore, complement protein results for stored samples need to be assessed against reference intervals determined under similar conditions.

Procedure

Notes

- Consult your BN[®] System Instruction Manual for details regarding operation of the instrument.
- Allow reagents and samples to equilibrate to room temperature (+15 to +25 °C) before use on a BN[®] 100 System. With a BN[®] II or BN ProSpec[®] System, reagents and samples stored at +2 to +8 °C can be used immediately.

Assay Protocols for BN[®] Systems

The assay protocols for serum and plasma are given in the BN[®] System Instruction Manual and software of the instrument. All steps are performed automatically by the system.

Establishment of the Reference Curve

Reference curves are generated by multi-point calibration. Serial dilutions of N Protein Standard SL are automatically prepared by the instrument using N Diluent. The standard dilutions are to be used within four hours.

The reference curves can be used for as long as controls with corresponding method-dependent target values, e.g. N/T Protein Control SL/L, M and H, are reproduced within their respective confidence intervals. If a different lot of antiserum is used, a new reference curve must be generated.

Assay of Specimens

Serum and plasma samples are automatically diluted 1:20 with N Diluent. The diluted samples must be measured within four hours. If the results obtained are outside the measuring range, the assay can be repeated using a higher or lower dilution of the sample. Refer to BN[®] System's Instruction Manual for information on repeat measurements using other dilutions.

Internal Quality Control

Assay N/T Protein Controls SL/L, M and H after each establishment of a reference curve, the first use of an antiserum vial as well as with each run of samples. The controls are assayed and evaluated as for patient samples. Assigned value and confidence interval are listed in the Table of Assigned Values of the respective control.

If the result of a control is outside the confidence interval, the determination must be repeated. If the repeated determination confirms the deviation, a new reference curve must be established. Do not release patient results until the cause of deviation has been identified and corrected.

Results

Evaluation is performed automatically in g/L or in a derived unit selected by the user on the BN[®] System.

Limitations of the Procedure

No interference was detected for concentrations of triglycerides up to 5.7 g/L (C3) or 2.4 g/L (C4), of bilirubin up to 600 mg/L, and of free hemoglobin up to 10 g/L. No interference from commonly used drugs is known. Turbidity and particles in the samples may interfere with the determination. Therefore, samples containing particles must be centrifuged prior to testing.

Dade Behring has validated use of these reagents on various analyzers to optimize product performance and meet product specifications. User defined modifications are not supported by Dade Behring as they may affect performance of the system and assay results. It is the responsibility of the user to validate modifications to these instructions or use of the reagents on analyzers other than those included in Dade Behring Application Sheets or these instructions for use.

Results of these tests should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.

Due to matrix effects, interlaboratory survey samples and control samples may yield results that differ from those obtained with other methods. It may therefore be necessary to assess these results in relation to method-specific target values.

Reference Intervals

The following reference intervals apply for serum and plasma samples from healthy adults⁵:

C3/C3c	0.9 - 1.8 g/L
C4/C4c	0.1 - 0.4 g/L

Fresh samples can be expected to have lower C3 concentrations. The reference intervals for C3 and C4 will depend on sample age and storage (refer to "Specimens" section).

Therefore, each facility should determine its own reference intervals since values may vary depending on the individual population studied and age and storage of the sample.

Specific Performance Characteristics

Assay Range and Analytical Sensitivity

The N Antisera to Human C3c and C4 are designed to measure complement factor concentrations within a range of approximately

C3c	0.12 - 4.1 g/L
C4	0.06 - 1.9 g/L

for a sample dilution of 1:20 in case of serum or plasma.

The analytical sensitivity of the assays is determined by the lower limit of the reference curve and depends therefore upon the concentrations of the proteins in the N Protein Standard SL. A typical sensitivity in serum or plasma for C3/C3c is 0.02 g/L and for C4/C4c 0.001 g/L.

Specificity

There are no known cross-reactivities of the antisera used.

Precision

The following coefficients of variation (CV) were obtained with N Antisera to Human Complement Factors on a BN[®] System (n = 40):

C3/C3c	mean (g/L)	run-to-run CV (%)	within-run CV (%)	total CV (%)
N/T Protein Control, Level L	0.84	2.0	2.4	2.9
N/T Protein Control, Level M	1.17	2.6	4.2	4.5
N/T Protein Control, Level H	1.65	2.8	2.8	3.8
Serum Pool (low)	1.21	2.1	3.2	3.5
Serum Pool (high)	1.71	2.0	2.0	2.7
C4/C4c	mean (g/L)	run-to-run CV (%)	within-run CV (%)	total CV (%)
N/T Protein Control, Level L	0.132	1.3	1.9	2.2
N/T Protein Control, Level M	0.197	1.9	1.9	2.6
N/T Protein Control, Level H	0.317	1.6	1.8	2.3
Serum Pool (low)	0.338	2.8	1.7	3.2
Serum Pool (high)	0.462	1.6	2.3	2.6

The results were generated by analysis of variance.

Method Comparison

Serum samples were assayed with N Antisera to Human Complement Factors on a BN[®] System (y) and an immunoturbidimetric method (x) Correlation of the results yielded the following data:

Protein	n	y (BN [®]) = 0.89 x - 0.07 g/L	Correlation Coefficient
C3/C3c	146	y (BN [®]) = 0.89 x - 0.07 g/L	0.93
C4/C4c	129	y (BN [®]) = 1.02 x - 0.04 g/L	0.95

Note

The values cited for specific performance characteristics represent typical values and are not to be regarded as specifications for the N Antiserum to Human C3c or C4.

Bibliography

1. Thomas L. The complement system. In: Thomas L, ed. Clinical Laboratory Diagnostics. Frankfurt: TH-Books Verlagsgesellschaft, 1998: 794-806.
2. Hebert LA, Cosio FG, Neff JC. Diagnostic significance of hypocomplementemia. Kidney Int 1991; 39: 811-21.
3. West CD. The complement profile in clinical medicine. Inherited and acquired conditions lowering the serum concentrations of complement component and control proteins. Complement Inflamm 1989; 6: 49-64.
4. Okumura N, Nomura M, Tada T, et al. Effects of sample storage on serum C3 assay by immunonephelometry. Clin Lab Sci 1990; 3: 54-57.
5. Dati F, Schumann G, Thomas L, et al. Consensus of a group of professional societies and diagnostic companies on guidelines for interim reference ranges for 14 proteins in serum based on the standardization against the IFCC/BCR/CAP Reference Material (CRM 470). International Federation of Clinical Chemistry. Community Bureau of Reference of the Commission of the European Communities. College of American Pathologists. Eur J Clin Chem Clin Biochem 1996; 34: 517-20.

BN ProSpec is a registered trademark of Dade Behring Marburg GmbH in the USA, Germany and other countries.

* BN is a trademark of Dade Behring Marburg GmbH in the USA.

Dade Behring Marburg GmbH
Emil-von-Behring-Str. 76
35041 Marburg/Germany
www.dadebehring.com

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