N Latex CDT Kit

OPCS G03 C0543 (1484)

Intended Use

In vitro diagnostic for the quantitative determination of carbohydrate-deficient transferrin (CDT) in human serum by means of particle-enhanced immunonephelometry using the BN[®] II and BN ProSpec[®] System for the identification of chronic elevated alcohol consumption.

Summary and Explanation

The main component of serum transferrin is a glycoprotein with a molecular mass of about 80 kD that contains two polysaccharide side chains. Each side chain carries two terminal sialic acid residues. Human serum transferrin can, however, occur in different isoforms with different degrees of glycosylation¹.

The tetrasialoform with two carbohydrate side chains, each featuring two terminal sialic acid residues, dominates in about 90% of healthy subjects¹. Transferrin molecules with only one carbohydrate side chain (disialotransferrin) or even no carbohydrate chain (asialotransferrin) occur in increased numbers as a result of alcohol consumption². These and other structurally altered transferrin isoforms induced by alcohol consumption are referred to as CDT (carbohydrate-deficient transferrin)^{1,3}.

As a rule, the daily alcohol consumption of about 50 – 60 g ethanol over two weeks leads to increased CDT levels. Increased CDT values can be expected to normalize after abstinence from alcohol for about two to four weeks, depending on the height of the CDT level^{3,4}.

The determination of CDT provides a valuable contribution to the recognition of patients with chronically high alcohol consumption, the monitoring of changes in alcohol consumption and the monitoring of abstinence. Different studies have shown that CDT is one of the specific markers for increased alcohol consumption over a longer period of time^{5,6}.

Non-alcohol induced illnesses that may trigger an increase in CDT include chronic active hepatitis, primary biliary cirrhosis, liver failure and the extremely rare CDG (carbohydrate-deficient glycoprotein) syndrome^{1,4}.

The calculation of %CDT from CDT and transferrin concentrations enables the influences of transferrin levels, iron status and mild to moderate limitations of liver function on the CDT result to be largely minimized. An additional improvement in diagnostic accuracy can be achieved by the combination of CDT and γ GT (γ -glutamyltransferase)^{5,7}.

Principle of the Method

Competitive assay

The CDT in the sample competes with CDT-coated polystyrene particles for the bond to specific monoclonal antibodies against human CDT, which are likewise bound to polystyrene particles. In the presence of CDT in the sample, there is no or little aggregation of the polystyrene particles. In the absence of CDT in the sample, the polystyrene particles aggregate. The higher the CDT content in the assay, the lower the scattered light signal. The evaluation is performed by comparison with a standard of known concentration.

Reagents

Materials provided

N Latex CDT Kit, Code No. OPCS

N CDT Reagent 1, three vials containing 0.9 mL each

N CDT Reagent 2, three vials containing 0.9 mL each

N CDT Supplementary Reagent, three vials containing 2 mL each

N CDT Standard SL, three vials containing 1 mL each

N CDT Control SL/1, three vials containing 1 mL each

N CDT Control SL/2, three vials containing 1 mL each

Composition and Standardization

N CDT Reagent 1 consists of a suspension of polystyrene particles coated with carbohydrate-deficient transferrin (< 0.017 g/L).

N CDT Reagent 2 consists of a suspension of polystyrene particles coated with monoclonal anti-CDT antibodies (mouse) (< 0.009 g/L).

N CDT Supplementary Reagent consists of a phosphate buffered saline solution with polyethylene glycol sorbitan monolaurate (5 g/L) and EDTA (< 0.2 mol/L).

N CDT Standard SL (human) consists of a stabilized human serum matrix and carbohydrate-deficient transferrin. The concentration of CDT was calibrated with reference to purified CDT. The analytical value for the N CDT Standard SL is indicated in the enclosed table.

N CDT Control SL/1 and N CDT Control SL/2 (human) each consist of a stabilized human serum matrix and carbohydrate-deficient transferrin. The concentration of CDT was calibrated with reference to the N CDT Standard SL. The assigned value and confidence interval for N CDT Control SL/1 and N CDT Control SL/2 are given on the enclosed table. These also contain the transferrin concentration and a %CDT value for the N CDT Controls SL/1 and SL/2; the transferrin concentration for these was calibrated with reference to N Protein Standard SL and the %CDT value was calculated from the respective CDT and transferrin concentration.

Preservative

N CDT Reagent 1 and N CDT Reagent 2: Gentamicin 6.25 mg/L, Amphotericin 0.625 mg/L

N CDT Supplementary Reagent, N CDT Standard SL, N CDT Control SL/1 and N CDT Control SL/2: Sodium azide < 1 g/L.

Warnings and Precautions

- For *in vitro* diagnostic use.
- Products containing sodium azide must be handled with due caution:
 - Do not ingest or allow to contact with skin or mucous membranes. When disposing of in water, flush with large amounts of water. Sodium azide can form explosive azides when contacting heavy metals such as copper or lead.
- Each individual blood donation intended for use in the manufacture of N CDT Reagent 1 and N CDT Reagent 2 as well as N CDT Standard SL, N CDT Control SL/1 and SL/2, was tested for HBsAg, Anti-HCV, Anti-HIV1 and Anti-HIV2. Only donations with negative findings were used for manufacture. Nevertheless, since the absence of infectious agents cannot be proven, all materials obtained from human tissues or body fluids should always be handled with due care, observing the precautions recommended for biohazardous material⁸.

Preparation of Reagents

The reagents, standard and controls are ready to use and require no additional preparation. Shake carefully to mix before the first use.

Storage and Stability

Store the reagents protected from light!

Stability at +2 to +8 °C:

The expiry date is indicated on the label;

Stability once opened:

Two weeks for all reagents, standard and controls if stored at +2 to +8 °C securely capped immediately after each use. Do not freeze the reagents.

Stability on-board the BN[®] II and BN ProSpec[®] Systems:

A minimum of three days, at eight-hours per days or a comparable period of time.

Note: "On-board" stability may vary, depending on the BN[®] System used and laboratory conditions. For further details, refer to BN[®] II and/or BN ProSpec[®] System Instruction Manual.

"On-board" stability of the N CDT Controls SL/1 and SL/2 on the BN ProSpec[®] System are stated in the Instruction Manual of the system.

Material required, but not provided

BN[®] II or BN ProSpec[®] System

N Supplementary Reagent L, Code No. OQTD

N Diluent, Code No. OUMT

BN[®] II Evaporation Stoppers (optional), Code No. OVLE

Optional, for the transferrin determination:

N Antiserum to Human Transferrin, Code No. OSAX

N Protein Standard SL, Code No. OQIM

N/T Protein Control SL/L, M, H, Code Nos. OQIN, OQIO, OQIP (for the quality assurance of the transferrin determination)

N Reaction Buffer, Code No. OUMS

Additional materials and supplies as described in your BN[®] System Instruction Manual.

Specimens

Suitable samples are human serum as either as fresh as possible (stored for no more than seven days at +2 to +8 °C) or stored frozen. Samples can be stored at below -20 °C for up to three months if they are frozen within 24 hours after collection and if repeated freeze-thaw cycles are avoided. Serum samples must be fully coagulated and must not contain any particles or traces of fibrin after centrifuging. Lipemic samples or frozen samples that became turbid after thawing must be clarified by centrifuging (10 minutes at approximately 15,000 x g) prior testing.

Procedure

Notes

- Consult your BN[®] System Instruction Manual for details regarding operation of the instrument.
- Reagents and samples stored at +2 to +8 °C may be used immediately for testing on a BN[®] II or a BN ProSpec[®] System.
- In order to receive the %CDT value on the BN[®] II and the BN ProSpec[®] System for a sample in addition to the CDT concentration the transferrin concentration of the respective sample must be determined with the transferrin assay available on the BN[®] Systems.
- The transferrin content given for each of the N CDT Controls SL/1 and SL/2 aids in the calculation of the %CDT of individual control measurements. It is not suitable for the quality control of the transferrin assay on the BN[®] Systems.

Assay Protocol on the BN[®] Systems

The assay protocol for serum is given in the BN[®] System Instruction Manual and software of the instrument. All steps are performed automatically by the system.

Establishment of the Reference Curve

Reference curves are generated by multi-point calibration. Serial dilutions of N CDT Standard SL are automatically prepared by the instrument using N Diluent. The standard dilutions must be used within four hours. The analytical value is indicated in the enclosed table. For the BN ProSpec[®] System, a lot data CD (Code No. OVLP) can be used to enter the analytical value. The reference curve is valid for two weeks. It can be used beyond this period of time as long as controls with corresponding method-dependent target values, e. g. N CDT Control, Level SL/1 and Level SL/2 are reproduced within their respective confidence interval. If a different lot of reagent is used, a new reference curve must be generated. The exact measuring range depends on the protein concentration of each lot of N CDT Standard SL. Typical measuring ranges are provided in the respective BN[®] Systems Instruction Manual.

Assay of Specimens

Samples are automatically diluted 1:5 with N Diluent. They must be measured within four hours. The manual selection of another sample dilution than 1:5 is not permissible.

Internal Quality Control

Assay N/T CDT Controls SL/1 and SL/2 after each establishment of a reference curve, the first use of a reagent vial, as well as with each series of samples. The controls are to be assayed and evaluated as for patient samples. The assigned value and confidence interval can be found in the enclosed table. For the BN ProSpec[®] System, a lot data CD (Code No. OVLP) can be used to enter the assigned value. The stated assigned values are intended for use as an intra-laboratory control for the assessment of precision and analytical bias. If used for precision control, the user should establish the target concentration and confidence limits during a preliminary phase. If a control value is outside the confidence interval, the determination must be repeated. If the repeated determination confirms the deviation, a new reference curve should be established. Do not release patient results until the cause of deviation has been identified and corrected.

Results

Evaluation is performed automatically in mg/L or in a unit selected by the user on the BN[®] System.

In addition, the calculation of %CDT has been integrated into the software of the BN[®] Systems. In the case of simultaneous determination of CDT and transferrin from the same sample, %CDT can be received additionally.

Limitations of the Procedure

Turbidity and particles in the sample may interfere with the determination. Therefore, samples containing particles must be centrifuged prior to testing. Lipemic samples or samples that contain particles that cannot be clarified by centrifuging (10 minutes at approx. 15,000 x g) must be excluded from testing.

Patient samples may contain heterophile antibodies that could react in immunoassays to give a falsely elevated or depressed result. This assay has been designed to minimize interference from heterophilic antibodies. Nevertheless, complete elimination of this interference from all patient specimens cannot be guaranteed.

Dade Behring has validated use of these reagents on various analyzers to optimize product performance and meet product specifications. User defined modifications are not supported by Dade Behring as they may affect performance of the system and assay results. It is the responsibility of the user to validate modifications to these instructions or use of the reagents on analyzers other than those included in Dade Behring Application Sheets or these instructions for use.

Results of these tests should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.

Due to matrix effects, inter-laboratory survey samples and control samples may yield results that differ from those obtained with other methods. It may therefore be necessary to assess these results in relation to method-specific target values.

Reference Interval

In a study of serum samples from healthy 561 adults from Central Europe in whom an elevated level of alcohol consumption had been excluded, the results determined with the N Latex CDT Reagent ranged from 28.1 - 76.0 mg/L CDT (1st to 99th percentile). The %CDT values in this population, with reference to the results obtained with the N Antiserum to Human Transferrin, ranged from 1.19 - 2.47 %CDT (1st to 99th percentile).

In addition, each laboratory should determine its own reference ranges since the results may vary depending on the population tested.

Specific Performance Characteristics

Sensitivity

The analytical sensitivity of the test is determined by the lower limit of the reference curve and is therefore dependent on the concentration of the analyte in the N CDT Standard SL. A typical detection limit for CDT is 20 mg/L.

Specificity

There are no known cross-reactions with the antibodies used.

Precision

The precision of N Latex CDT was calculated by measuring the controls and two serum pools with different CDT concentrations and calculating the coefficients of variation (CV).

Sample	Precision (n = 40)			
	Mean value	Within the series	From run to run	Total
	CV	CV	CV	CV
N CDT Control SL/1	61.6	4.2	1.6	4.2
N CDT Control SL/2	166.2	2.8	1.5	3.0
Serum Pool 1	49.4	4.9	7.6	8.9
Serum Pool 2	213.9	3.5	2.4	4.0

The coefficient of variation for %CDT of these samples was calculated based on the result obtained with the N Antiserum to Human Transferrin using the BN[®] System and ranged from 1.8% to 9.8%.

Method comparison

The %CDT values determined for serum samples from 200 test subjects were measured, on the one hand, with the aid of high-performance liquid chromatography (HPLC), on the other hand using a BN[®] System. If a cut-off of 2.5 %CDT is set for %CDT values obtained on the BN[®] System, the result is a specificity of 97 % and a sensitivity of 93 % in comparison to the %CDT values determined with HPLC.

Note

The values cited for specific performance characteristics represent typical results and are not to be regarded as specifications for the N Latex CDT Kit.

Bibliography

- Golka K, Wiese A. Carbohydrate-deficient transferrin (CDT) - a biomarker for long-term alcohol consumption. J Toxicol Environ Health B 2004;7:319-37.
- Peter J, Unverzagt C, Engel WD, et al. Identification of carbohydrate deficient transferrin forms by MALDI-TOF mass spectrometry and lectin ELISA. Biochim Biophys Acta 1998; 1380: 93-101.
- Helander A. Biological markers in alcoholism. J Neural Transm 2003; 66 (Suppl): 15-32.
- Allen JP. Use of biomarkers of heavy drinking in health care practice. Mil Med 2003; 168: 364-7.
- Anttila P, Järvi K, Latvala J, et al. A new modified γ -CDT method improves the detection of problem-drinking: studies in alcoholics with or without liver disease. Clin Chim Acta 2003; 338: 45-51.
- Schwan R, Albuissou E, Malet L, et al. The use of biological laboratory markers in the diagnosis of alcohol misuse: an evidence-based approach. Drug Alcohol Depend 2004; 74: 273-9.
- Chen J, Conigrave KM, Macaskill P, et al. Combining carbohydrate-deficient transferrin and gamma-glutamyltransferase to increase diagnostic accuracy for problem drinking. Alcohol Alcoholism 2003; 38: 574-82.
- U.S. Department of Health and Human Services CDC, Biosafety in Microbiological and Biomedical Laboratories, HHS Publication (CDC) 93-8395; 1999; Section II; 8-16.

* BN is a trademark of Dade Behring Marburg GmbH in the USA.

BN ProSpec is a registered trademark of Dade Behring Marburg GmbH in the USA, Germany and other countries.

Dade Behring Marburg GmbH
Emil-von-Behring-Str. 76
35041 Marburg/Germany
www.dadebehring.com

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