Complement Reagents

for the kinetic determination of the activity of the total complement

Intended Use

The complement reagents permit investigation of the functional capacity of the entire complement system. The complement components covered are those of the classical pathway as well as of the alternative pathway. Inhibitory components of the complement system are not covered. Used in combination with the appropriate deficient plasma or serum, the reagents permit quantitative measurement of the activity of individual components of the complement system.

Summary and Explanation

The complement system is an endogenous plasmatic defence system. It is activated either by antigen-antibody complexes (classical pathway) or by foreign surfaces (alternative pathway). Pathological conditions may result in a reactive increase in the concentrations of the individual components (acute-phase reaction). Diminished activities are found as a result of consumption or congenital deficiencies.

Like the CH50 test, the method measures the functional activity of the complement. The test is therefore suitable for detecting hereditary defects as well as acquired deficiencies (e.g. due to immune complex diseases such as glomerulonephritis, systemic lupus erythematosus, generalized vasculitis or cryoglobulinaemia) and for monitoring their course.

Principle of the Method

Incubation of the sample with sensitized erythrocytes leads to activation of the classical complement pathway via C1q. C3b formed by C3 convertase becomes a component of C5 convertase as well as an initiator of the amplification reaction with the involvement of Factor B and Factor D. This process also results in the formation of alternative C5 convertase. C5b released by C5 initiates the lysis reaction with the involvement of the remaining Factors C6 to C9.

The lysis of the erythrocytes leads to a decrease in the optical density. The time taken for the lysis of a defined quantity of erythrocytes is used as a basis for determining the complement activity. Deficiencies or diminished levels of individual components (e.g. due to prior activation of the complement system) result in a prolongation of the time taken for lysis. Shorter times are obtained when the levels of complement proteins are increased in reaction to corresponding causative events.

Reagents

Materials provided

Complement Reagents, Test kit, Code No. OWZD with

5 x for 5 mL Erythrocyte Reagent
1 x 25 mL Amboceptor Reagent

Composition:

Erythrocyte Reagent: Sheep erythrocytes in stabilizer solution (sodium chloride approximately 9 g/L, HEPES 2.4 g/L).

Preservatives: 5-Chloro-2-methyl-4-isothiazole-3-one (0.37 g/L),
2-methyl-4-isothiazole-3-one (0.13 g/L).

Amboceptor Reagent: Rabbit antibodies to sheep erythrocytes, sodium chloride, calcium chloride and magnesium chloride, HEPES (2.4 g/L).

Preservative: Sodium azide (0.5 g/L).

Warnings and Precautions

- For *in vitro* diagnostic use.
- As complete absence of infectious agents can never be assured, all samples (e.g. patient plasma) and products (e.g. standard and control preparations) derived from human blood should be handled with due care, observing the precautions recommended for biohazardous material³.
- Reagents containing sodium azide must be handled with due caution: Do not ingest or allow to contact skin or mucous membranes! Sodium azide can form explosive azides when contacting heavy metals such as copper or lead.

Preparation of Reagents

Before the first use, carefully resuspend the erythrocytes by means of gentle swirling. Watch for erythrocytes adhering to the stopper! Next mix the Erythrocyte Reagent with the labelled volume of Amboceptor Reagent and leave to stand at +37 °C for approx. 45 min. Shake briefly before use and bring to +37 °C.

Storage and Stability

Stored unopened at +2 to +8 °C the test kit can be used by the expiry date given on the label.

Stability of the freshly prepared reagent mixture:

+37 °C 4 hours
+2 to +8 °C 24 hours

Do not freeze the reagent!

The Erythrocyte Reagent must exhibit a bright red colour. If the reagent has developed a brown colour it must not be used in the assay.

Different coagulation analyzers have individual stability data.

Materials required but not provided

Standard Human Plasma, Code No. ORKL

Control Plasma N, Code No. ORKE

Control Plasma P, Code No. OUPZ

Coagulation analyzer (see Chapter titled "Limitations of the Procedure")

Equipment

Complement Reagents can be used in a great number of coagulation analyzers.

Observe the operating instructions issued by the analyzer's manufacturer!

Specimens

Citrated plasma or serum.

Stability of the sample:

-20 °C 1 month
+15 to +25 °C 4 hours

Samples stored at -20 °C are to be thawed at +37 °C within 10 min and then assayed within 2 hours.

Procedure

Manual:

Cuvette: 1 cm
Wave length: 578 nm (500 - 600 nm)
Test temperature: +37 °C
Recorder: 1 to 5 cm/minute

Prewarm plastic cuvettes to +37 °C. Used the sample in undiluted form.

Sample Reagent	50 µL 1000 µL
Determination of the lysis time: Measure the time taken for the absorbance at 578 nm to decrease by 0.100.	

Procedure for the Chromotimer: Use half the volumes of sample and reagent. The Chromotimer automatically measures the lysis time and calculates the final result.

Automatic:

Specific instructions for coagulation analyzers may be obtained from Dade Behring by request.

Evaluation

The results are evaluated using a reference curve prepared by serial dilution of a fresh plasma pool or standard plasma, e.g. Standard Human Plasma, with isotonic saline:

Complement activity	Standard plasma µL	Isotonic saline µL
100 %	100	0
75 %	75	25
50 %	50	50
25 %	25	75
10 %	10	90

In order to calibrate for 150 or 200 and 300 % the volume of the standard can be reduced to 75, 100 and 150 µL respectively. The lysis times of the standard series are determined in the above assay and plotted against the activity in percent of the norm on log-log graph paper.

The complement activity of the sample is read from the reference curve. If the labelled assigned value of the Standard Human Plasma is not 100 % of the norm, but e.g. only 95 %, the result read is multiplied by 0.95.

In order to improve the precision in the range below an activity of 50 %, the volume of sample used can be doubled. In this case the sample must be evaluated via a calibration curve prepared using the same ratio of sample to reagent.

Serum samples can also be evaluated from a calibration curve prepared using plasma.

The calibration curve must be re-determined for each change of device and for each new lot of Complement Reagents.

Internal Quality Control

A Control must be measured with each calibration and in each series of measurements. Suitable controls are Control Plasma N and Control Plasma P. The controls should be processed just like the samples. Each laboratory should determine its own quality control range, either by means of the target values and ranges provided by the manufacturer of the controls or by means of its own confidence level established in the laboratory. If the measured control value lies outside the confidence level previously established, then the reagents, the calibration curve and the coagulation analyzer should be examined.

Limitations of the Procedure

There are appropriate coagulation analyzers for which Dade Behring provides operating instructions. The utilization of Complement Reagents in other coagulation analyzers must be validated by the testing laboratory. In doing so, performance data other than that provided may be obtained.

EDTA-containing samples must not be used in the assay.

Reference Interval

70 - 140 % of the norm

Systematic deviations from these ranges may be due to the instrument used. If necessary, the laboratory will have to establish its own reference intervals.

Specific Performance Characteristics

Measuring Range

The measuring range extends from 10 to 300 % of the norm.

Specificity

The test covers the classical as well as the alternative pathways of the complement system but not the inhibitory components.

Precision

For normal plasma samples the intra-assay coefficients of variation range between 2.7 and 6.5 % for manual pipetting and from 0.8 to 2.5 % for automatic pipetting. The inter-assay CV was found to be 4.6 (manual method) and 1.7 % (automatic method).

For pathological plasma samples the corresponding intra-assay values ranged between 4.3 and 14 % or 0.8 and 5.3 % respectively. The inter-assay CV was found to be 5.5 (manual method) and 2.4 % (automatic method).

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Dade Behring Marburg GmbH